

التاريخ: 14/12/2025

Q079.2026

رقم الاستدراج:

الشركة:

تلفون:

إيميل:

BR2526-152

استدراج عروض أسعار: الرجاء تزويدنا بأسعار اللوازم التالية شاملة لضريبة القيمة المضافة:

اسم الصنف	الوحدة	الكمية	السعر	اجمالي تكلفة البند	مدة التوريد
Disintegration Tester	عدد	1			
Friability Tester	عدد	1			
Tablet Hardness Tester	عدد	1			
Dissolution Tester	عدد	1			
المجموع الكلي (ش.ض)			شيقل		

ختم وتوقيع المورد:

تعليمات الاستدراج

** آخر موعد لتسليم العروض هو يوم:	الثلاثاء
** الموافق:	23/12/2025
** طبيعة البنود المطلوبة:	أجهزة
** بعملة:	شيقل
** الأسعار تشمل التوصيل ل:	الحرم الرئيسي - جنين
** نوع التوريد:	توريد كامل
** يرجى إرسال العرض ب:	بالظرف المغلق

شروط الاستدراج:

** ضرورة تحديد مدة الكفالة إن وجدت.
** ضرورة تحديد مدة التوريد من تاريخ استلام امر الشراء، وفي حال عدم الالتزام يتم فرض غرامة تأخير .
** ختم وتوقيع المورد المعتمد على هذا النموذج، وكتابة رقم الاستدراج واسم الشركة كما هو مسجل لدى الضريبة على الظرف.
** للمراجعة والاستفسار بخصوص المواصفات يرجى الاتصال على هاتف 042418888 رقم داخلي 1488.
** للمراجعة والاستفسار بخصوص المواصفات عن طريق فاكس 042510972 أو بريد الكتروني samah.zaidalkelani@aaup.edu
** لمشاهدة العطاءات والاستدراجات، يرجى الدخول الى الرابط http://www.aaup.edu/tenders



Item #1: Disintegration Tester

Brand name.

3-year warranty, 10-year availability of spare parts.

Origin: Brand name

- tablet disintegration tester.
 - Compliant with USP <701> and Ph. Eur. 2.9.1 for disintegration testing.
 - With two baskets, a clear view, and easy cleaning.
 - With magnetic coupling.
 - Quick lock
 - Digital heart (independent).
 - Low liquid level.
 - Basket for large tablets and capsules
 - Clear display
 - Complete system with all accessories and ready to use
 - السعر يشمل التوريد والتركيب ويجب أن يتم توضيح جميع متطلبات التركيب في عرض السعر من ملحقات ومتطلبات.
- ❖ *Kindly your offer will not be taken into our consideration if you don't fill in the table below.*

Disintegration Tester	Yes/No	Comment
Manufacture	Manufacture	
Model	Origin	
Origin:	Model	
Delivery time	Time	
3 years Full warranty,10 years availability for spare part	Yes / No	
tablet disintegration tester.	Yes / No	
Compliant with USP <701> and Ph. Eur. 2.9.1 for disintegration testing.	Yes / No	
With two baskets, a clear view, and easy cleaning.	Yes / No	
With magnetic coupling.	Yes / No	
Quick lock	Yes / No	
Digital heart (independent).	Yes / No	
Low liquid level.	Yes / No	
Basket for large tablets and capsules	Yes / No	
Clear display	Yes / No	
Complete system with all accessories and ready to use	Yes / No	
Price includes installation training on-site	Yes / No	
User's and service manual	Yes / No	
Date of manufacture of the machine not less than 2025	Yes / No	



Item #2 Friability Tester:

Brand name.

3 years warranty, 10 years availability of spare parts.

Origin: USA, EURO, Japan

- Tablet friability tester for pharmaceutical laboratories.
- Compliant with USP <1216> and Ph. Eur. 2.9.7 for friability testing.
- Drum made of transparent material for easy observation of tablets.
- Digital timer with adjustable rotation speed.
- Capable of holding standard tablet weights.
- Quick and easy drum mounting and removal.
- Low noise operation.
- Clear digital display for speed, time, and test parameters.
- Complete system with all accessories, delivered ready to use.
- Operation and service manuals included.
- Installation and training included.

السعر يشمل التوريد والتركيب ويجب أن يتم توضيح جميع متطلبات التركيب في عرض السعر من ملحقات ومتطلبات.

❖ **Kindly your offer will not be taken into our consideration if you don't fill in the table below.**

Friability Tester	Yes/No	Comment
Manufacture.	Manufacture	
Model.	Model	
Origin:	Origin	
Delivery time.	Weeks	
3 years Full warranty.	Yes/No	
Brand name.	Yes/No	
Tablet friability tester for pharmaceutical laboratories.	Yes/No	
Compliant with USP <1216> and Ph. Eur. 2.9.7 for friability testing.	Yes/No	
Drum made of transparent material for easy observation of tablets.	Yes/No	
Digital timer with adjustable rotation speed.	Yes/No	
Capable of holding standard tablet weights.	Yes/No	
Quick and easy drum mounting and removal.	Yes/No	
Low noise operation.	Yes/No	
Clear digital display for speed, time, and test parameters.	Yes/No	
Complete system with all accessories, delivered ready to use.	Yes/No	
Operation and service manuals included.	Yes/No	
Installation and training included.	Yes/No	
Complete and Ready to use with all accessories.	Yes/No	
10 years availability for spare parts.	Yes/No	
Date of manufacture not less than 2025.	Yes/No	
في حال وجود أية قطعة ضرورية لتشغيل الجهاز، غير موضحة يرجى ذكرها وتوضيحها في العرض غير ذلك تعتبر الشركة بتحمل وجود أية نواقص.	Yes/No	



Item #3 Tablet Hardness Tester**Brand name.****3 years warranty, 10 years availability of spare parts.****Origin: USA, EURO, Japan**

- Tablet hardness tester for pharmaceutical laboratories.
 - Compliant with **USP <1217>** and **Ph. Eur. 2.9.8** for hardness testing.
 - Suitable for tablets of different shapes and sizes.
 - Automatic and manual operation modes.
 - Digital display for hardness, force units (N, kg/cm², or lb), and test parameters.
 - Quick and easy calibration system.
 - Data storage and export capability (optional).
 - Low noise operation.
 - Complete system with all accessories, delivered ready to use.
 - Operation and service manuals included.
 - Installation and training included.
 - ومتطلبات. السعر يشمل التوريد والتركيب ويجب أن يتم توضيح جميع متطلبات التركيب في عرض السعر من ملحقات
- ❖ **Kindly your offer will not be taken into our consideration if you don't fill below table.**

Tablet Hardness Tester	Yes/No	Comment
Manufacture.	Manufactur e	
Model.	Model	
Origin:	Origin	
Delivery time.	Weeks	
3 years Full warranty.	Yes/No	
Brand name.	Yes/No	
<input type="checkbox"/> Tablet hardness tester for pharmaceutical laboratories.	Yes/No	
<input type="checkbox"/> Compliant with USP <1217> and Ph. Eur. 2.9.8 for hardness testing.	Yes/No	
<input type="checkbox"/> Suitable for tablets of different shapes and sizes.	Yes/No	
<input type="checkbox"/> Automatic and manual operation modes.	Yes/No	
<input type="checkbox"/> Digital display for hardness, force units (N, kg/cm ² , or lb), and test parameters.	Yes/No	
<input type="checkbox"/> Quick and easy calibration system.	Yes/No	
<input type="checkbox"/> Data storage and export capability (optional).	Yes/No	
<input type="checkbox"/> Low noise operation.	Yes/No	
<input type="checkbox"/> Complete system with all accessories, delivered ready to use.	Yes/No	
<input type="checkbox"/> Operation and service manuals included.	Yes/No	
User, Application, and service manuals.	Yes/No	
Price includes installation with application training on site.	Yes/No	
Complete and Ready to use with all accessories.	Yes/No	
Date of manufacture not less than 2025.	Yes/No	
في حال وجود أية قطعة ضرورية لتشغيل الجهاز غير موضحة يرجى ذكرها وتوضيحها في العرض غير ذلك تعتبر الشركة بتحمل وجود أية نواقص.	Yes/No	



Item #4: Dissolution Tester**Brand name.****3-year warranty, 10-year availability of spare parts.****Origin: Brand name**

- Tablet/capsule dissolution tester for pharmaceutical laboratories.
 - Compliant with USP <711> and Ph. Eur. 2.9.3 for dissolution testing.
 - Six or more vessels with precise temperature control.
 - Transparent vessel for clear observation.
 - Digital display for temperature, rotation speed (RPM), and test duration.
 - Magnetic or mechanical coupling for paddle/basket rotation.
 - Automatic detection of low liquid level and overheating.
 - Quick vessel and paddle/basket mounting/removal.
 - Stirring speed adjustable and programmable.
 - Data output and recording capability (optional).
 - Complete system with all standard accessories, delivered ready to use.
 - Operation and service manuals included.
 - Installation and training included.
- السعر يشمل التوريد والتركيب ويجب أن يتم توضيح جميع متطلبات التركيب في عرض السعر من ملحقات ومتطلبات.
- ❖ **Kindly your offer will not be taken into our consideration if you don't fill in the table below.**

Dissolution Tester	Yes/No	Comment
Manufacture	Manufacture	
Model	Origin	
Origin:	Model	
Delivery time	Time	
3 years Full warranty,10 years availability for spare part	Yes / No	
Tablet/capsule dissolution tester for pharmaceutical laboratories.	Yes / No	
Compliant with USP <711> and Ph. Eur. 2.9.3 for dissolution testing.	Yes / No	
Six or more vessels with precise temperature control.	Yes / No	
Transparent vessel for clear observation.	Yes / No	
Digital display for temperature, rotation speed (RPM), and test duration.	Yes / No	
Magnetic or mechanical coupling for paddle/basket rotation.	Yes / No	
Automatic detection of low liquid level and overheating.	Yes / No	
Quick vessel and paddle/basket mounting/removal.	Yes / No	
Stirring speed adjustable and programmable.	Yes / No	
Data output and recording capability (optional).	Yes / No	
Complete system with all standard accessories, delivered ready to use.	Yes / No	
Price includes installation training on-site	Yes / No	
User's and service manual	Yes / No	
Date of manufacture of the machine not less than 2025	Yes / No	

